



Polaris Medical Inc.

www.polarismedicalinc.com



SHOWER TROLLEY



GENIUS

Shower trolley – The easiest way to shower persons who need assistance.



- Easy care through operator-friendly construction.
- Sturdy and maintenance-free.
- Running and steering wheels, diam. 125 mm, lockable.
- Lying area properly dimensioned for patients, practically designed components.
- Skin-friendly shower cover and neck cushion with thermally insulated plastic.
- Side rail can be easily folded down and locked.
- Insertable, flexible drain hose.
- Comfortable adjustment to the desired working height through maintenance-free foot pump which can be operated from either side.
- Couch can be hinged up for easy cleaning.
- Shower cover easily removed though large Velcro fastener.
- Couch of stainless steel, surface can be sterilised.
- Chassis powder-coated, anti-slip and impact-resistant.

GENIUS

Distributed by:
Polaris Medical Group

T 905.467.5152
F 905.333.0407
TF 877-377-9988

info@polarismedicalinc.com
www.polarismedicalinc.com



005 Prog. No. 154119
EN ISO 9001